



## Nasal Washer

Model: NW1



Patent Nos:  
ZL 2014 3 0313271.3  
Other patents pending

### • It's Effective

Flushes nasal cavities and rhinopharynx to clear secretions, mucus and harmful substances and collects them in a chamber.

### • It's Efficient

Nasal washer nebulizes 10 ml of solution in about 2 minutes.

### • It's Comprehensive

Nasal Washer is applicable to the whole range of Rossmax compressor nebulizers.

### • It's Comfortable

Silicone nasal adaptors provide comfortable and soft touch and reliable seal.

### • It's Convenient

Easy to use with saline solutions and warm waters or medications commonly available in pharmacies.

### • It's Non-Toxic

Made of biocompatible, non-allergenic materials and Latex and BPA free.


### • It's Practical

Comes with 3 sizes of nasal adaptors and suitable for the whole range.

### Technical Specification

- Fluid chamber capacity: Max 20ml
- Single patient use only
- Applicable models: NB500, NA100, NE100, NB80, NB60, NI60, NK1000, NJ100

\* Subject to technical modification without prior notice.



Model	Qty per carton	Carton volume
NW1	50 pcs	0.031 cbm/ctn
QTY per 20'	QTY per 40'	QTY per 40'HQ
45,200 pcs	90,350 pcs	104,850 pcs

● Storage and Transportation Condition:  
-20°C~60°C(14°F~140°F)  
10%~90% RH, 700~1060 hPa